

EW80[®] *mat*

- **Aldehyde-free**
- **Chlorine-free**
- **Alcohol-free**
- **pH-neutral**



CE 0482

Compatibility with materials checked:



EW80[®]
mat

Cleaning and disinfection of PPE with EW80[®] *mat*

EW80[®] *mat* has been tested by leading manufacturers and is particularly suitable for the mechanical and thermochemical treatment of respiratory protective equipment and chemical protection suits made of water- and temperature-compatible materials.

The short application time of the 1% EW80[®] *mat* solution ensures that the equipment is reliably treated in just five minutes. Due to low water consumption per disinfection cycle, the environment is sustainably relieved.

EW80[®] *mat* is tested in accordance with the VAH/DGHR guidelines for instrument disinfection and, in line with EU regulations as per the Medical Devices Act, CE marked.

Art.-no.

510205
510230

packaging

3 x 5 Liter canister
30 Liter canister